

2007: PRO Order Packet

Application Form and Professional Use Waiver



MTS Roller™

MTS Roller™ (Microneedle Therapy System™) is a “meso-therapy” skin treatment device designed and manufactured by **Clinical Resolution, Inc.** (a California corporation). All MTS roller models are currently available for U.S. market (see below) and are FDA registered as Type I supplementary medical instruments. All components are produced in high-tech ISO9001 quality-controlled manufacturing facilities located in Seoul, South Korea and Switzerland. MTS Roller is patented.

General Description

MTS Roller is a durable and high-quality precision tool. The barrel and handle are made from Lexan® polycarbonate resin. The axle is polished stainless steel axle. The 200 bristles (needles) – the highest count of any roller – are made from surgical stainless steel. A patented manufacturing process holds the needles permanently in place to within 0.02mm of the specified length from the roller base. All rollers have unique quality-control serial numbers.

GUARANTEE: 100% Contaminant-Free

MTS Rollers are sterilized with gamma radiation and sealed in a hermetic pouch and are guaranteed to be 100% contaminate-free upon delivery.

Storage and Use

Polycarbonate resin does not degrade from prolonged use or storage. Personal use rollers (CR2, CR3) should be replaced after about 6 months of regular use (2 to 3 times week). Clinical rollers (CR5, CR10, CR20, CR30) should be replaced after 10 applications. Rollers that have been dropped or bumped should be replaced immediately. Per FDA rules, rollers should not be transferred between patients. Keep rollers in the protective case when not in use. For more information, consult CRL protocol guides.

Care of MTS Rollers™

Current MTS Roller models can withstand repeated immersion in “cold sterilization” agents such as germicides (phenol, ethyl alcohol) and Cidex™ (glutaraldehyde). Rollers should be completely dry before soaking to prevent dilution of the sterilization agent. Soak time depends upon the manufacturer's specifications. Rinse rollers thoroughly before using and return to its protective case after drying.



Pictured: Model CR30 (3.0mm)

MTS ROLLER™ MODELS:

- CR2:** 0.2mm (personal use, “Stage 1”)
- CR3:** 0.3mm (personal use, “Stage 2”)
- CR5:** 0.5mm (clinical use only)
- CR10:** 1.0mm (clinical use only)
- CR20:** 2.0mm (clinical use only)
- CR30:** 3.0mm (clinical use only)

CLINICAL RESOLUTION

3255 Wilshire Blvd., Suite 1708
Los Angeles, CA 90010

Copyright © 2007 by Clinical Resolution Lab, Inc.

2007: PRO Order Packet

Application Form and Professional Use Waiver



To order MTS Roller™

Complete the forms contained in this "Pro Package" as a convenient means to begin a unique opportunity for professional clinicians (physicians and/or other registered/licensed/certified skincare professionals) to test MTS Roller in their practice at DEALER PRICES (see "Pro Price List"). PLEASE FAX BACK ALL COMPLETED FORMS TO OUR OFFICE - **213-365-0053** – or email a PDF version of the completed document to your sales representative. Please note that there are TWO portions to this Packet: "Certificate of Professional Clinical Status," and "Waiver of Liability". Call your sales representative if you have questions.

Certification of Professional Clinical Status

Practitioner's Name (Buyer):

License/Certification Type:

Are you currently in practice? YES NO

Business Name (dba):

Mailing/Billing Address (stamp or label):

Your Clinic Specialty(ies)

Shipping Address (if different than mailing/billing):

Years in practice

Reseller Permit # (if in California)

CRL ADMIN:

Phone:

Fax:

Email/Website

I am the "Practitioner" named above. All of the information I have provided is accurate to the best of my ability. (Sign and date)

CLINICAL RESOLUTION 

3255 Wilshire Blvd., Suite 1708
Los Angeles, CA 90010

Copyright © 2007 by Clinical Resolution Lab, Inc.

2007: PRO Order Packet

Application Form and Professional Use Waiver



Waiver of Liability

Buyer (named signor on the "Certificate of Professional Clinical Status" and also signing below) assumes sole responsibility for the use of MTS rollers. Buyer also assumes complete and full responsibility – without exception – for any liability (material, financial, medical, and/or psychological) that may arise from the use of MTS rollers and completely indemnifies and holds harmless Clinical Resolution Lab, Inc., its officers, representatives, and/or agents from any judgment, any arbitration, and/or any other settlement brought with or without cause.

Clinical Resolution Lab, Inc. will replace defective MTS rollers that fail during normal usage. Buyer also acknowledges that until further notice, Clinical Resolution advises that MTS Rollers should be sterilized by cold sterilization processes between usages. **Do not autoclave MTS rollers.** This warranty is therefore exempt for MTS Rollers that have been autoclaved.

MTS Rollers CR5, CR10, CR20 and CR30 are sold as single-patient use supplemental surgical tools for use in the clinical environment only under supervision of a physician. MTS Rollers CR2 and CR3 are sold as single-patient use take-home tools to be used under the guidance of a physician.

Other terms and conditions apply per CRL, Inc. invoicing and payment policies.

I agree to the terms and conditions of the above Waiver of Liability (sign and date)



3255 Wilshire Blvd., Suite 1708
Los Angeles, CA 90010

Copyright © 2007 by Clinical Resolution Lab, Inc.